

A Paradigm Shift in Hemodynamic Monitoring: The Expanding Value of Venous Oximetry

INTRODUCTION

Once reserved for only the most critically ill patients, hemodynamic monitoring has evolved to play a significant role in the management of a wider range of patient populations.¹ Advanced technologies provide comprehensive data, including venous oxygen saturation, to guide proactive management of patients prior to hemodynamic crisis rather than reacting to late indicators of instability.

While venous oximetry, specifically central venous oximetry (ScvO₂), has most commonly been associated with Early Goal-Directed Therapy (EGDT) for the identification and treatment of severe sepsis and septic shock, it has also been clearly linked with improving outcomes in patients experiencing decreased cardiac output due to myocardial infarction, cardiac surgery, trauma or hemorrhagic shock, high-risk surgeries, and respiratory failure.^{2,3,4,5,6,7}

This paper outlines the utility of venous oximetry for a variety of medical conditions as established through medical literature, clinical benefits and challenges of implementing it, and an overview of monitoring technologies that can be considered.

WHY MONITOR?

The purpose of hemodynamic monitoring is to ensure adequate tissue perfusion and oxygenation. This process occurs naturally in a healthy body as it responds to increased demand for oxygen or reduced arterial oxygen by increasing cardiac output to establish homeostasis. In the presence of critical illness or injury, however, the heart's ability to respond may be limited. As a result, depleted tissue draws oxygen from the venous oxygen reserve, leading to anaerobic metabolism and the production of lactic acid. If that is not addressed promptly, lactic acidosis and global tissue hypoxia ensues.

In recent years, improvements in survival rates for acute myocardial infarction, trauma, stroke, and sepsis have been realized through early identification of tissue oxygen imbalance using hemodynamic monitoring combined with tissue perfusion therapy. Additional evidence indicates that early initiation of hemodynamic resuscitation consistently improves mortality rates.⁸

WHAT TO MONITOR?

Continuous Measurement of SvO₂ & ScvO₂

The principles of oxygen transport and utilization are simple. Tissue consumes a percentage of delivered oxygen from blood for cellular respiration. The blood then returns to the venous circulation, where its oxygen content can be measured with a blood sample from the pulmonary artery or from the superior vena cava. Oxygenation measurements taken at the pulmonary artery are called mixed venous oxygen saturation, or SvO₂, while measurements taken at the superior vena cava are called central venous oxygen saturation, or ScvO₂. When a patient is unable to generate cardiac output sufficient to meet the metabolic needs of the tissue, the tissue extracts greater amounts of oxygen, leaving less oxygen present in the venous blood. This is reflected as a decreasing SvO₂ or ScvO₂ level. Conversely, when tissue does not extract the oxygen, both SvO₂ and ScvO₂ levels rise.

Improvements in survival rates for acute myocardial infarction, trauma, stroke, and sepsis have been realized through early identification of tissue oxygen imbalance using hemodynamic monitoring combined with tissue perfusion therapy.

Various studies have shown that there can be slight differences between ScvO₂ and SvO₂ values—with either measurement reading 5-6% higher or lower than the other based on patient acuity.^{9,10,11,12} However, the two have been shown to maintain a high correlation (see Figure 1).⁹ As a result, ScvO₂ is accepted as an appropriate surrogate for SvO₂.

In post-cardiac arrest patients, ScvO₂ may detect ongoing oxygen imbalances that do not fully resolve during resuscitation. ScvO₂ values below 40-50% indicate the likelihood of re-arrest, and ScvO₂ values above 60-70% indicate hemodynamic stability.²⁰

WHO BENEFITS?


Continuous venous oximetry can play a vital role in the effective management of critical care patients. Beyond identifying early signs of severe sepsis, ScvO₂ monitoring is an effective tool for diagnosing and managing patients during high-risk surgeries and conditions of myocardial infarction, cardiac failure, respiratory failure, and trauma (see Table 1).^{2,3,4,5,13}

Sepsis

Sepsis is the condition most closely associated with use of continuous venous oximetry as part of EGDT protocols and the most studied to-date. In

a study comparing continuous ScvO₂ monitoring with intermittent monitoring as a means to guide use of EGDT in sepsis, continuous monitoring had a greater success rate in achieving an ideal ScvO₂ target of 70% within six hours of diagnosis, resulting in quicker improvement of global tissue perfusion and improved survival.¹⁴ Another study of continuous ScvO₂ monitoring in the treatment of sepsis found a 16% reduction in in-hospital mortality, a reduction in length of stay of approximately four days, and a savings of \$12,000 per discharge.¹⁵

FIGURE 1. TIME COURSE OF MIXED AND CENTRAL VENOUS O₂ SATURATION DURING DIFFERENT EXPERIMENTAL PERTURBATIONS. HTS=HYPERTONIC SALINE SOLUTION (7.5%)⁹



A recent study evaluating EGDT using ScvO₂ has shown that the EGDT group and use of ScvO₂ resulted in outcomes comparable to usual-care in hospitals experienced in managing the septic population. Though the EGDT-managed group's total cost was higher when using a CVC oximetry catheter for ScvO₂, the authors of the study stated the cost difference was not significant.¹⁶ Regardless, EGDT protocols have had an influential impact on clinical outcomes in the last decade, and the use of continuous venous oximetry monitoring is an important part of that protocol.¹⁷

In addition, patients at risk for severe sepsis often present with comorbidities where the availability of venous oximetry provides clinicians with invaluable physiological data with which to more effectively manage their patients.

Cardiac Care

Nearly five million Americans live with congestive heart failure (CHF), resulting in roughly 875,000 hospitalizations each year where the primary diagnosis is CHF.¹⁸ Although cardiac patients are at high risk for developing global tissue hypoxia, vital signs and findings from physical examination have been the primary diagnostic tools guiding emergent care. Research supporting the use of hemodynamic monitoring via ScvO₂ in this patient population is growing. The use of continuous ScvO₂ monitoring, coupled with lactic acid measurement, has been shown to be superior to the assessment of vital signs alone for identifying tissue hypoxia in cases of cardiogenic shock and for guiding the implementation of early, aggressive management strategies (see Table 1).³

Several studies have also shown the usefulness of ScvO₂ monitoring for patients in acute cardiac arrest.^{6,19} One study reported that with an ScvO₂ level greater than 60%, return of spontaneous circulation was likely to occur, while an ScvO₂ value greater than 72% was associated with a certain return of spontaneous circulation.³ Additional research suggests that in post-cardiac arrest patients, ScvO₂ may detect ongoing oxygen imbalances that do not fully resolve during resuscitation. Specifically, an ScvO₂ value below 40-50% indicates the likelihood of re-arrest, and an ScvO₂ value above 60-70% indicates that the patient is hemodynamically stable.²⁰

Roughly 10% of cardiac patients require prolonged postoperative care because of hemodynamic instability, organ dysfunction, or multiple organ failure—and this patient population is likely to grow as surgical and technological advances enable surgery on traditionally high-risk patients.²¹ In the immediate post-cardiopulmonary bypass period, inadequate oxygen delivery is a known predictor of a prolonged ICU stay, while the recognition of improved cardiac output and subsequent enhancement of oxygen consumption are clearly associated with improved outcomes.^{22,23} The information gleaned from ScvO₂ monitoring also helps to guide the delivery of intravenous fluid and vasopressor/inotropic therapy for these cardiac patients.²⁴

TABLE 1. CLINICAL APPLICATIONS AND BENEFITS OF VENOUS OXIMETRY MONITORING (SvO₂ AND ScvO₂)^{2,3,4,5,13}

PATIENT TYPE	
> Septic	> Myocardial Infarction (MI)
> Post-Coronary Artery Bypass Graft Surgery (CABG)	> Trauma
> High-Risk Surgery	> Respiratory Failure
BENEFIT OF MONITORING	
> Indication of adequate cardiac output (CO)	> Indication of O ₂ consumption
> Decreased length of stay (LOS)	> Indication of additional interventional needs
> Decreased morbidity/mortality	> Improved outcomes

For patients with trauma and hemorrhagic shock, an ScvO₂ value less than 65% measured in the first 24 hours after admission is associated with higher mortality and prolonged ICU and hospital length of stay.²⁹

Trauma

Traumatic injury may result in tissue hypoxia, oxidative stress, capillary leakage, and metabolic changes.²⁵ More importantly, tissue hypoxia within the first 24 hours after a traumatic event has been shown to be a predictor of multiple organ dysfunction syndrome.²⁶ Yet, as with many other critical care situations, the assessment of vital signs remains the gold standard used to detect and treat shock states in trauma patients, when in fact, ScvO₂ monitoring has demonstrated significant advantages over vital signs assessment.

For patients with trauma and hemorrhagic shock, an ScvO₂ value less than 65% may indicate that additional resuscitation or surgical intervention is needed.^{27,28} In addition, an ScvO₂ value less than 65% measured in the first 24 hours after admission is associated with higher mortality and prolonged ICU and hospital length of stay.²⁹ ScvO₂ measurement was also shown to be a trustworthy parameter to estimate blood loss, especially in patients with an ScvO₂ level below 65% despite stable vital signs.⁴

This lack of linear correlation between low ScvO₂ and stable vital signs has been corroborated by a non-randomized study of patients presenting to the emergency department in shock conditions who presented an ScvO₂ value below 65% in almost half of cases, despite other hemodynamic parameters appearing stable.⁵ Conversely, patients with ScvO₂ values greater than 65% discharged alive from the ICU had significantly shorter ICU and total hospital lengths of stays.²⁹


High-Risk Surgery

Surgical patients are subject to the effects of anesthesia, surgical trauma, inflammation, fluid shifts, and blood loss. In extensive surgical procedures these effects are magnified, which may impair oxygen delivery to organs and tissues, triggering organ failure. The development of postoperative organ failure in 27-77% of high-risk surgeries, such as organ

transplants, severely effects the prognosis of patients and substantially increases the utilization of resources and cost of care.³⁰ Length of stay and postoperative mortality are also increased in these high-risk patients.^{7,31} Consequently, the use of early and efficient means to detect and treat potential triggers of organ failures, such as tissue hypoxia, is particularly important in this population.

In a study that used ScvO₂ and arterial lactate as two indices of tissue level oxygenation, lactate consistently rose later than ScvO₂.³⁰ Interestingly, organ failures were observed much more often in patients with at least one elevated lactate value. Thus, the researchers concluded that staving off organ failure is reliant on earlier and more aggressive hemodynamic management before lactic acidosis commences.

FIGURE 2. VALUES OF SvO₂ DURING MECHANICAL VENTILATION AND A TRIAL OF SPONTANEOUS BREATHING IN THE SUCCESS GROUP³⁸



The rapid decrease of SvO₂ at the onset of the spontaneous breathing trial is highly predictive of unsuccessful weaning.


Respiratory failure

There are as many as 200,000 cases of acute respiratory distress syndrome (ARDS) and acute lung injury (ALI) reported in the U.S. each year.³² Due to the increased risk of complications, the use of pulmonary arterial catheterization for routine monitoring of these patients is not recommended.³³ However, the World Health Organization advises healthcare providers to “closely monitor patients with severe acute respiratory infections for signs of clinical deterioration, such as severe respiratory distress/respiratory failure or tissue hypoperfusion/shock, and apply supportive care interventions.”³⁴ In the early phase of ARDS, an associated septic state is usually responsible for fluid retention, and hemodynamic optimization is critical to stabilizing the patient.³⁵ Therefore, ScvO₂ is especially important in managing this complicated patient population.

ARDS patients require longer periods of mechanical ventilation which are associated with longer ICU stays.³⁶ Weaning ARDS patients from the ventilator can be a difficult and prolonged process. Generally, tolerance of a spontaneous breathing test (SBT) indicates weaning success, but less effectively predicts extubation success. In fact, the need for reintubation within 24 to 72 hours occurs in 5-30% of patients, dependent on the population.³⁷ Research shows that patients who failed weaning also failed to increase oxygen delivery to tissues (see Figure 2).³⁸ This respiratory exhaustion, a common challenge post-extubation, is reflected in near real-time through ScvO₂ measurement. As such, in the weaning process, measurement of ScvO₂ has the potential to be a reliable and convenient tool to rapidly warn of acute changes in the oxygen supply and demand of these patients.³⁹ In separate studies, researchers reported that a reduction in ScvO₂ by $\geq 4.5\%$ during the SBT was an independent predictor of reintubation.⁴⁰

Another study showed the difference between ScvO₂ at the 1st and 30th minute of a spontaneous breathing trial as a predictor of extubation failure. Their multivariate regression identified this difference (Δ ScvO₂) as the only independent variable able to discriminate extubation outcomes. The receiver operating characteristics curve yielded an area under the curve (AUC) of 0.861 (see Figure 3).⁴¹

FIGURE 3. ROC CURVE OF Δ ScvO₂ IN PREDICTING EXTUBATION FAILURE (AUC=0.861, 95% CI 0.702 - 0.954, P < 0.001)⁴¹



The AUC=0.861 shows that Δ ScvO₂ has a very good discriminative predictability for when a patient will fail to extubate.

In the early phase of ARDS, an associated septic state is usually responsible for fluid retention, and hemodynamic optimization is critical to stabilizing the patient.³⁵ Therefore, ScvO₂ is especially important in managing this complicated patient population.

VENOUS OXIMETRY MONITORING TECHNOLOGIES

Pulmonary Arterial Catheters (PAC)

SvO₂ is not readily measurable without a pulmonary artery catheter, a technology that has created enthusiasm as well as controversy in its 40-year history. Apart from its pressure and SvO₂ monitoring capabilities, pulmonary arterial catheterization can measure cardiac output using a method called thermodilution. Cardiac output measurements from PACs are then combined with other hemodynamic measurements to calculate systemic and pulmonary vascular resistance. Use of PACs does not come without risk of complications, including pneumothorax, hemorrhage, pulmonary artery rupture, and pulmonary embolism,⁴² leading some clinicians to question their value.⁴³

Central Venous Oximetry Catheters (CVOC)

ScvO₂ measurements are traditionally obtained by a CVC with fiber optic oximetry, as opposed to a more invasive PAC. Placement of a CVC is easier to perform and has fewer inherent risks than placement of a PAC, allowing ScvO₂ monitoring to benefit patients and clinicians in a variety of settings where insertion of a PAC is not feasible.¹¹

There has been considerable discussion regarding ScvO₂ as a substitute for SvO₂, and an extensive bibliography supports its use in hemodynamic monitoring.^{44,45} At the center of this research, a prospective cohort study found the venous-arterial oxygenation differences obtained from the pulmonary artery and central venous circulations were equal and inversely correlated with the cardiac index.⁴⁵

Another study showed that an ScvO₂ value may be used as a surrogate marker to reflect adequacy or inadequacy of cardiac output.⁴⁶ Studies suggest ScvO₂ may show that a lower cardiac output or mean arterial pressure (MAP) is sufficient to meet the current metabolic demands of the patient, reducing reliance on pressors, oxygen/mechanical ventilation, blood products and/or inotropes.^{47, 48}

Utilizing a PICC line for venous oximetry would allow for less invasive and more cost-effective ScvO₂ monitoring without the infection potential associated with a direct neck puncture.

While less invasive than PAC placement, the use of CVCs is also not without risk. Complications of CVC insertion may include pneumothorax, hematoma, carotid artery puncture, and infection. Additionally, because a CVC is intended for a short period of use, repeated insertions are necessary when long-term monitoring is required, impacting patient comfort and increasing the risk of infection while potentially limiting future vascular access.

Less Invasive Venous Oximetry

Peripherally inserted central venous catheters (PICCs) have been safely used for many years to access the central venous circulation to administer fluids such as parenteral nutrition, chemotherapy, vasopressor, antibiotics and other solutions, and facilitate blood draws. More than one million PICCs are placed each year in the U.S. alone.⁴⁹ Because of this ability to access the central venous circulation at the superior vena cava, PICCs are a less invasive avenue to consider for venous oximetry.

Currently, the only available PICC line with oximetry monitoring technology is the TriOx[®]-PICC minimally invasive oximetry sensor (ICU Medical Inc., San Clemente CA). TriOx-PICC uses three-wavelength oximetry technology to filter noise and artifact caused by cell orientation, vessel wall reflections, and changes in pH.

The use of a PICC has several significant advantages over the use of CVCs. PICCs avoid the complications associated with the direct neck puncture of the central venous circulation system.^{50,51,52} These devices have also been shown to have lower catheter-related bloodstream infection (CRBSI) rates as compared to CVCs⁵³ and support a longer period of use.⁵⁴

Because PICC lines can be inserted by trained nurses outside the operating room, PICC use reduces cost and decreases treatment delays common in busy operating or intensive care units.⁵⁵ So in addition to being less invasive and potentially safer, PICCs are also a more convenient, efficient, and cost-effective method to access ScvO₂ levels.

CONCLUSION

Venous oximetry has been established as a valuable method to monitor and maintain adequate tissue oxygenation and avoid hemodynamic crisis in critically ill and at-risk patients. ScvO₂ has been shown to be a surrogate measurement for SvO₂, reducing the reliance on invasive pulmonary artery catheterization. Its broad utility has been demonstrated across a variety of conditions, from cardiac care to trauma, respiratory failure and high-risk surgery, in addition to its origins in the treatment of sepsis.

The technology used to obtain venous oximetry measurements has evolved from pulmonary artery catheterization to less invasive central venous catheterization, with the goal of increasing patient safety, reducing costs, and minimizing procedure time. As the industry continues to migrate toward less invasive and more cost-effective patient monitoring trends, PICC line oximetry will continue to become a viable and efficient method for monitoring ScvO₂.



The TriOx-PICC minimally invasive oximetry sensor (ICU Medical, Inc., San Clemente CA) allows clinicians to monitor continuous ScvO₂ less invasively and more cost-effectively than with traditional central lines.

References

- Hamilton MA, Cecconi M, Rhodes A. A Systematic Review and Meta-Analysis on the Use of Preemptive Hemodynamic Intervention to Improve Postoperative Outcomes in Moderate and High-Risk Surgical Patients. *Anesth Analg* 2011;112:1392–402.
- Astiz ME, Rackow EC, Kaufman B, Falk JL, Weil MH. Relationship of oxygen delivery and mixed venous oxygenation to lactic acidosis in patients with sepsis and acute myocardial infarction. *Crit Care Med* 1988, 16:655–658.
- Ander DS, Jaggi M, Rivers E, Rady MY, Levine TB, Levine AB, Masura J, Gryzbowski M. Undetected cardiogenic shock in patients with congestive heart failure presenting to the emergency department. *Am J Cardiol* 1998, 82:888–891.
- Scalea TM, Hartnett RW, Duncan AO, Atweh NA, Phillips TF, Sclafani SJ, Fuortes M, Shaftan GW. Central venous oxygen saturation: a useful clinical tool in trauma patients. *J Trauma* 1990, 30:1539–1543.
- Rady MY, Rivers EP, Martin GB, Smithline H, Appelton T, Nowak RM. Continuous central venous oximetry and shock index in the emergency department: use in the evaluation of clinical shock. *Am J Emerg Med* 1992, 10:538–541.
- Nakazawa K, Hikawa Y, Saitoh Y, et al. Usefulness of central venous oxygen saturation monitoring during cardiopulmonary resuscitation: a comparative case study with end-tidal carbon dioxide monitoring. *Intensive Care Med.* 1994; 20:450–451.
- Shoemaker WC, Appel PL, Kram HB. Role of oxygen debt in the development of organ failure sepsis, and death in high risk surgical patients. *Chest* 1992; 102:208–215
- Otero RM, Nguyen HB, Huang DT, et al. Early goal-directed therapy in severe sepsis and septic shock revisited: concepts, controversies, and contemporary findings. *Chest* 2006; 130(5):1579–95.
- Reinhart K, Rudolph T, Bredle D, et al. Comparison of Central-Venous to Mixed-Venous Oxygen Saturation During Changes in Oxygen Supply/Demand. *Chest.* 1989 Jun; 95(6):1216–1221.
- Shahbazi S, Khademi S, Shafa M, et al. Serum Lactate Is Not Correlated with Mixed or Central Venous Oxygen Saturation for Detecting Tissue Hypo Perfusion During Coronary Artery Bypass Graft Surgery: A Prospective Observational Study. *Int Cardiovasc Res J.* 2013;7(4):130–134.
- Reinhart K, Huhn HJ, Hartog C, et al. Continuous central venous and pulmonary artery oxygen saturation monitoring in the critically ill. *Intensive Care Med.* 2004;30:1572–1578.
- Squara P. Central venous oxygenation: when physiology explains apparent discrepancies. *Critical Care* 2014, 18:579.
- Krafft P, Steltzer H, Hiesmayr M, Klimscha W, Hammerle AF. Mixed venous oxygen saturation in critically ill septic shock patients. The role of defined events. *Chest* 1993, 103:900–906.
- Ising P, Smith TW and Simpson SQ. Effect of Intermittent vs. Continuous ScvO₂ Monitoring On Sepsis Bundle Compliance And Mortality. October 2009, Vol 136, No. 4. Meeting Abstracts
- Rivers E, Nguyen B, Havstad S, et al. Early goal-directed therapy in the treatment of severe sepsis and septic shock. *N Engl J Med.* 2001;345:1368–1377
- Mouncey PR, Osborn TM. Trial of Early, Goal-Directed Resuscitation for Septic Shock. *N Engl J Med.* 2015 Apr 2;372(14):1301–11
- Levy MM, Dellinger RP, Townsend SR, et al. The Surviving Sepsis Campaign: results of an international guideline-based performance improvement program targeting severe sepsis. *Critical Care Medicine* 2010;38:367–74.
- Go AS et al. AHA Statistical Update. Heart Disease and Stroke Statistics—2014 Update. *Circulation* 2014; 129: e28–e292
- Snyder AB, Salloum LJ, Barone JE, et al. Predicting short-term outcome of cardiopulmonary resuscitation using central venous oxygen tension measurements. *Crit Care Med.* Jan 1991; 19:111–113.
- Rivers EP, Ander DS, Powell, D. Central venous oxygen saturation monitoring in the critically ill patient. *Curr Opin Crit Care.* 2001; 7:204–211.
- Higgins TL, Yared JP, Ryan T. Immediate postoperative care of cardiac surgical patients. *J Cardiothorac Vasc Anesth* 1996;10:643–58
- Routsis C, Vincent JL, Bakker J, De Backer D, Lejeune P, d'Hollander A, et al. Relation between oxygen consumption and oxygen delivery in patients after cardiac surgery. *Anesth Analg* 1993; 77:1104–10.
- Polonen P, Ruokonen E, Hippelainen M, Poyhonen M, Takala J. A prospective, randomized study of goal-oriented hemodynamic therapy in cardiac surgical patients. *Anesth Analg* 2000; 90:1052–9.
- Shoemaker WC, Appel PL, Waxman K, Schwartz S, Chang P. Clinical trial of survivors' cardiorespiratory patterns as therapeutic goals in critically ill postoperative patients. *Crit Care Med* 1982; 10:398–403.
- Keel M, Trentz O. Pathophysiology of polytrauma. *Injury* 2005, 36:691–709.
- Blow O, Magliore L, Claridge JA, Butler K, Young JS. The golden hour and the silver day: detection and correction of occult hypoperfusion within 24 hours improves outcome from major trauma. *J Trauma* 1999, 47:964–969.
- Wo CC, Shoemaker WC, Appel PL, et al. Unreliability of blood pressure and heart rate to evaluate cardiac output in emergency resuscitation and critical illness. *Crit Care Med.* 1993; 21:218–223.
- Kowalenko T, Ander D, Hitchcock R, et al. Continuous central venous oxygen saturation monitoring during the resuscitation of suspected hemorrhagic shock [abstract]. *Acad Emerg Med.* 1994; 1:A69.
- Di Filippo A, Gonnelli C, Perretta L, et al. Low central venous saturation predicts poor outcome in patients with brain injury after major trauma: a prospective observational study. *Scandinavian Journal of Trauma, Resuscitation and Emergency Medicine* 2009, 17:23
- Donati A, Loggi S, Preiser JC, et al. Goal-directed intraoperative therapy reduces morbidity and length of hospital stay in high-risk surgical patients. *Chest.* 2007 Dec; 132(6):1817–24.
- Lobo SM, Salgado PF, Castillo VG, et al. Effects of maximizing oxygen delivery on morbidity and mortality in high-risk surgical patients. *Crit Care Med* 2000; 28:3396–3404
- Rubinfeld GD, Caldwell E, Peabody E, et al. Incidence and outcomes of acute lung injury. *N Engl J Med* 2005; 353:1685–1693.
- Dellinger et al. Surviving Sepsis Campaign: International Guidelines for Management of Severe Sepsis and Septic Shock: 2012. *Critical Care Medicine* 2013; 40(2).
- World Health Organization. Clinical management of severe acute respiratory infections when novel coronavirus is suspected: What to do and what not to do. 2013. http://www.who.int/csr/disease/coronavirus_infections/InterimGuidance_ClinicalManagement_NovelCoronavirus_11Feb13u.pdf
- Stapleton RD, Wang BM, Hudson LD, et al. Causes and timing of death in patients with ARDS. *Chest* 2005, 128:525–532
- Ranieri VM, Rubinfeld GD, Thompson BT, Ferguson ND, Caldwell E, Fan E, Camporota L, Slutsky AS. Acute respiratory distress syndrome: the Berlin definition. *JAMA.* 2012; 307:2526–2533.
- Epstein SK. Decision to extubate. *Intensive Care Med.* 2002 May; 28(5):535–46.
- Jubran A, Mathru M, Dries D, Tobin MJ. Continuous recordings of mixed venous oxygen saturation during weaning from mechanical ventilation and the ramifications thereof. *Am J Respir Crit Care Med.* 1998; 158(6):1763–9.
- Sandham JD, Hull RD, Brant RF, Knox L, Pineo GF, Doig CJ, et al. A randomized, controlled trial of the use of pulmonary-artery catheters in high risk surgical patients. *N Engl J Med.* 2003; 348(1):5–14.
- Teixeira C, da Silva NB, Savi A, Vieira SR, Nasi LA, Friedman G, et al. Central venous saturation is a predictor of reintubation in difficult-to-wean patients. *Crit Care Med.* 2010 Feb; 38(2):491–6.
- Helmy TA, Baess AI, Abdelazez KMI. Central Venous Oxygen Saturation as a predictor of extubation failure in mechanically ventilated Chronic Obstructive Pulmonary Disease patients. *Journal of American Science* 2014;10(6)
- Hadian M and Pinsky MR. Evidence-based review of the use of the pulmonary artery catheter: impact data and complications. *Crit Care.* 2006; 10(Suppl 3):S8.
- Boldt J. Clinical review: hemodynamic monitoring in the intensive care unit. *Crit Care* 2002; 6(1):52–9.
- Dueck MH, Klimek M, Appenrodt S, Weigand C, Boerner U. Trends but not individual values of central venous oxygen saturation agree with mixed venous oxygen saturation during varying hemodynamic conditions. *Anesthesiology* 2005; 103: 249–57.
- Rigg JD, Nightingale PN, Faragher EB. Influence of cardiac output on the correlation between mixed venous and central venous oxygen saturation. *British Journal of Anaesthesia* 1993; 71: 459.
- Cuschieri J, Rivers EP, Donnino MW, Katilius M, Jacobsen G, Nguyen HB, et al. Central venous-arterial carbon dioxide difference as an indicator of cardiac index. *Intensive Care Med* 2005;31:818–22.
- Takala J. Hypoxemia due to increased venous admixture: influence of cardiac output on oxygenation. *Intensive Care Medicine.* 2007; 33(5):908–11.
- Walley KR. Use of central venous oxygen saturation to guide therapy. *American Journal of Respiratory and Critical Care Medicine.* 2011; 184(5):514–20.
- Moureau N. Vascular safety: it's all about PICCs. *Nurs Manage.* 2006; 37:22–27. Quiz 50.
- Maki DG, Kluger DM, Crnich CJ. The risk of bloodstream infection in adults with different intravascular devices: A systematic review of 200 published prospective studies. *Mayo Clin Proc* 2006; 81: 1159 – 71.
- O'Grady NP, Alexander M, Dellinger EP, et al. Guidelines for the prevention of intravascular catheter-related infections. *MMWR Recomm Rep.* 2002; 51(RR-10):1–29.
- Skies DJ, Abbott M, Keiser P. Peripherally inserted central catheters in patients with AIDS are associated with a low infection rate. *Clin Infect Diseases.* 2000; 30:949–952.
- Carrico R, ed. APIC Text of Infection Control and Epidemiology, 2nd ed. Washington, DC: Association for Professionals in Infection Control and Epidemiology; 2005.
- Johansson E, Hammarskjöld F, Lundberg D, Arnlind MH. Advantages and disadvantages of peripherally inserted central venous catheters (PICC) compared to other central venous lines: a systematic review of the literature. *Acta Oncol.* 2013 Jun; 52(5):886–92.
- Oakley C, Wright E, Ream E. The experiences of patients and nurses with a nurse-led peripherally inserted central venous catheter line service. *Eur J Oncol Nurs* 2000; 4: 207 – 18.